1D111

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EN Instruction for use

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NO Brukerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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	MULTI-AXIS ANKLE – 1D111 <i>User manual for orthopaedic prosthetists</i> Read before use	1D11199 2021-07
---	--	--------------------

Provide the patient with the instructions in § 3, 7, 8, 9

1. COMPONENTS INCLUDED

Designation	Reference	Included / Sold separately
Multi-axis ankle	1D111	Included
Front stop - Soft	EE010-1	Included
Front stop - Medium	EE010-2	Sold separately
Front stop - Hard	EE010-3	Sold separately
Rear shock absorber - Hard	EE009	Sold separately
Rear shock absorber - Medium	EE008	Included
Rear shock absorber - Soft	EE007	Sold separately
Rear shock absorber - Extra Soft	EE013	Sold separately
Setting kit	1D11150	Sold separately EE010-2, EE010-3 & EE009 included

2. DESCRIPTION, PROPERTIES AND ACTION MECHANISM

A. Description

Multi-axis ankle delivered assembled, and consisting of:

- Ankle body
- Front stop - Soft
- Rear shock absorber - Medium

B. Properties

Weight	250 g
Total height	65 mm
Design height	53 mm
Maximum patient weight (Including load borne)	100 kg

This device has been tested as per standard ISO 10328 for a load level of P5 (i.e. 100 kg) for 3 million cycles, corresponding to a service life of 4 to 5 years, depending on the patient's activity level.


C. Action mechanism

The multi-axis ankle adds the dorsal bending, plantar bending, inversion and eversion movements to any compatible foot.

3. INTENDED USERS/INDICATIONS

This medical device is supplied to health professionals (orthopaedic prosthetists), who will train the patient in their use. The prescription is drawn up by a doctor, alongside the orthopaedic prosthetist, to assess whether the patient is suited to using it.

 This device is for **SINGLE-PATIENT** use. It must not be reused on another patient.

This device is exclusively intended to be used as orthopaedic equipment for an above-knee or below-knee amputee with low-medium or normal activity level (2 or 3).

Maximum weight (including load borne): 100 kg

4. CLINICAL BENEFITS

The device will facilitate flattening the foot, and is suited to rough terrain, thanks to its multiple movement options.

A wide range of front stops and rear shock absorbers enables the ankle to be adapted to the patient's activity level and weight.


5. ACCESSORIES AND COMPATIBILITIES

Upper joint	Male pyramid
Lower joint	M10 screw

6. MOUNTING AND FITTING ON THE PATIENT

A. Assembly

- ⚠ Adhere to the tightening torque of the M10 screw: 22Nm
- ⚠ Adhere to the engaged screw thread length: 10 +/-1 mm

7. DETECTING MALFUNCTIONS

- ⚠ If you observe abnormal behaviour or feel modifications in the device's properties, or if it undergoes a heavy impact, please consult your orthopaedic prosthetist.

8. WARNINGS, CONTRAINDICATIONS, SIDE EFFECTS

A. Warnings

- ⚠ Using the device in a way that goes against the advice given by your orthopaedic prosthetist may cause deterioration of its components (e.g. bearing heavy loads, excessive strain, using beyond the service life, etc.).

- ⚠ The device is weatherproof only.

B. Contraindications

- ⚠ This device is not intended to be used for activities associated with a risk of significant impact or excessive loading.
- ⚠ This device is not resistant to fresh water, salt water or chlorinated water.

C. Side effects

There are no side effects directly linked to the device.

Any serious incident arising in relation to the device must be reported to the manufacturer and to the competent authority of the Member State.

9. CARE, STORAGE, DISPOSAL AND SERVICE LIFE


A. Care/ cleaning

Clean the device with a damp cloth or sponge.

- ⚠ No greasing, adjustments to the fasteners or any other operation is required for this device, except by the orthopaedic prosthetist.
- ⚠ After accidental contact with water:

- Rinse the device in clean water
- Dry thoroughly

The front stop and rear shock absorber may require replacement.

B. Storage

Usage and storage temperature: -20°C to +60°C

Relative air humidity: no restrictions

C. Disposal

The various components of the device are special waste: Aluminium, steel, thermoplastic, polyurethane elastomer, natural rubber. They must be treated in accordance with current, applicable legislation.

D. Service life

It is advisable to arrange for an orthopaedic prosthetist to perform an annual check.

10. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Risk identified		CE marking and year of 1 st declaration
---	--------------	---	-----------------	---	--

11. REGULATORY INFORMATION

This product is a CE marked medical device, and certified compliant with Regulation (EU) 2017/745

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PROTEOR SAS