
1M115

1M149

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EN Instruction for use

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NO Brukerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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
	BRAKE KNEE – 1M115 / 1M149 <i>User manual for orthopaedic prosthetists</i> Read before use	1M11599 2021-07
---	--	--------------------

Provide the patient with the instructions in § 3, 7, 8, 9

1. COMPONENTS INCLUDED

Designation	Reference	Included / Sold separately
Knee	1M115 or 1M149	Included
Extension return kit	M11510	Included and sold separately
Single-piece trim	1G13	Sold separately

2. DESCRIPTION, PROPERTIES AND ACTION MECHANISM

A. Description


Single-axis knee with stance phase brake and extension return.

Comes in two versions, and has 3 settings:

- Extension return force
- Joint friction
- Load-dependent brake sensitivity

B. Properties

Reference	1M115	1M149
Material	Stainless steel	Titanium
Weight	485 g	376 g
Height (H)	22 mm	30 mm
Total height (TH)	75 mm	
Design height (DH)	50 mm	
Maximum bend	150°	
Maximum patient weight (Including load borne)	100 kg	

This device has been tested as per standard NF EN ISO 10328 for a load level of P5 (i.e. 100 kg) for 3 million cycles, corresponding to a service life of 4 to 5 years, depending on the patient's activity level.

C. Action mechanism

This knee with stance phase brake comes with 3 distinct settings and extension return.

During the stance phase, the ground reaction force closes a jaw, which locks the rotation axis.

3. INTENDED USERS/INDICATIONS

This medical device is supplied to health professionals (orthopaedic prosthetists), who will train the patient in their use. The prescription is drawn up by a doctor, alongside the orthopaedic prosthetist, to assess whether the patient is suited to using it.

 This device is for SINGLE-PATIENT use. It must not be reused on another patient.

This device is designed only for the prosthetic system for a transfemoral amputee, or amputee with hip/knee disarticulation. It is specifically recommended for medium to low-activity patients (L1/L2).

Maximum weight (including load borne): 100 kg

 **Not suitable for children.**

 The knee's maximum bend is 150°. It may however be limited by the socket volume or by the aesthetic coating.

4. CLINICAL BENEFITS

The device makes it possible to:

- Set the brake according to the patient's weight and activity.
- Set the knee friction.
- Set the extension return.

5. ACCESSORIES AND COMPATIBILITY


Reference	1M115 / 1M149	
Upper joint	Tie-plate	1K176
	Connectors	1K172, 1K185, 1K66, 1K207
Lower joint	Tube and connector	1D41

6. MOUNTING AND FITTING ON THE PATIENT

A. Mounting the extension return

Install the extension return components: cable (a), settable stop (R), mobile support (c) and springs (d), as shown in the diagram.

- ⚠ If when the knee is bent, the cable (a) rubs against the side screws (2.1) (2.2) of the connector (2), replace these screws with the shortest screws (length 12 mm) contained in the package.


B. Alignments

Static alignments:

- ⚠ The performance of the braking mechanism depends primarily on the alignment of the prosthesis. In the sagittal plane, the load line runs via the Trochanter, between 5 and 15 mm in front of the knee axis.
- ⚠ The load line must run in front of the rear axis, or else it could cause premature deterioration of the knee, and lead to a risk of breakage and therefore falling.
- ⚠ Make sure to respect the patient's flexum. For the foot, comply with the alignment instructions from the manufacturer. In the frontal plane, the load line will run via the middle of the knee and the foot.

Dynamic alignments:

In normal walking, after fitting the extension return and the various settings set out below, the knee must be extended during the stance phase.

Check the alignment using an appropriate tool (laser, plumb line, etc.)

C. Setting

For the patient's safety, the orthopaedic prosthetist must perform the initial test with the factory settings, between parallel bars.

- ⚠ The foot model used may influence the knee settings. The knee should be set whenever the foot has been replaced. Using energy return feet is recommended. Ensure the patient's safety under all circumstances, particularly when they are going down slopes.
- ⚠ After any maintenance operation, the knee will need to be set again.

Then adjust the settings, preferably in the **order below**:

1. Setting the extension return force (R)

Set the extension return force by screwing in or unscrewing the settable stop (R) using a 4 mm male Allen key (k):

- *Clockwise = increase the tension*
- *Anticlockwise = decrease the tension.*

2. Setting the joint friction value (F)

The friction value has been factory set, but can be adapted according to the patient's approach. Use a 3 mm Allen key to make this setting, by adjusting screw (F):

- Clockwise = increase the friction
- Anticlockwise = decrease the friction


Rotating the key by 15 / 20° will make a perceptible change.

⚠ Too free a setting may lead to noise and premature wear. However, with too tight a setting, braking may occur inadvertently, even if the prosthesis is not under load.

3. Setting the load-dependent brake sensitivity (S)

The brake sensitivity can be adjusted via screw (S), using a 5 mm Allen key:

- Clockwise = joint brakes at higher load
- Anticlockwise = joint brakes at lower load.



D. Finishing

It is recommended to reinforce the trim before the joint, to extend its service life.

⚠ The trim can modify the initial settings of the knee, and limit certain movements. You should always check the operation of the knee once the trim has been mounted.

Thread the protective sleeve (g) onto the tube (h), and slide it upward until the stud (c1) on the mobile support (c) is inserted into the hole (g1) of the protective sleeve (g).

The socket must be in contact with the tube holding collar(3) at maximum bend. Any other contact zone is prohibited.

7. DETECTING MALFUNCTIONS

⚠ If you observe abnormal behaviour or feel modifications in the device's properties, or if it undergoes a heavy impact, please consult your orthopaedic prosthetist.

8. WARNINGS, CONTRAINDICATIONS, SIDE EFFECTS

A. Warnings

⚠ To avoid degrading the knee, do not use talc to eliminate friction noise, but a silicone spray instead. Talc causes deterioration of the mechanical components, which can lead to a malfunction with a risk of the patient falling.

PROTEOR disclaims all liability if the knee is used with talc.

⚠ There is a risk of the knee joint trapping a finger or catching clothing. To prevent any risk of injury due to joint movement, make sure that no-one puts their fingers near or inside the mechanism.

⚠ **The knee can withstand bad weather, but needs to be dried after being soaked.**

B. Contraindications

⚠ It is strictly prohibited to tighten or loosen any screw on this knee, with the exception of the four setting screws R, F and S, which are intended for use by the orthopaedic prosthetist.

⚠ Never grease the knee axes, as this could quickly degrade them.

- ⚠ The warranty excludes deterioration caused by misuse, unsuitable alignment, use in a very dusty environment and without suitable protection, or any inappropriate use.
- ⚠ Avoid exposing the knee to environments that could cause corrosion of the metal parts (soft water, seawater, chlorinated water, acids, etc.).
- ⚠ **It is prohibited to shower or bathe with the prosthesis**, as this could degrade its resistance and performance.
- ⚠ Never use this device near a heat source: risk of burns.
- ⚠ Use of solvents is prohibited.

C. Side effects

There are no side effects directly linked to the device.

Any serious incident arising in relation to the device must be reported to the manufacturer and to the competent authority of the Member State.

9. CARE, STORAGE, DISPOSAL AND SERVICE LIFE

A. Care/ cleaning

- ⚠ You can clean the knee using a wet sponge
- ⚠ Do not immerse it or expose it to water
- ⚠ After a storm (rain) or inadvertent splashing, the knee must be dried.

The extension return may require replacement, and can be ordered separately.

B. Storage

- ⚠ Usage and storage temperature: -10°C to +40°C
- Relative air humidity: no restrictions

C. Disposal

The various components of this device are special waste: elastomer, titanium, steel and brass. They must be treated in accordance with the locally applicable legislation.

D. Service life

It is advisable to arrange for an orthopaedic prosthetist to perform an annual check.

10. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Risk identified		CE marking and year of 1st declaration
---	--------------	---	-----------------	---	--

11. REGULATORY INFORMATION

This product is a CE marked medical device, and certified compliant with Regulation (EU) 2017/745

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PROTEOR SAS