# HYDEAL II

**1P120 / 1P120-KD**

**1P121 / 1P121-KD**

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**EN** Instruction for use

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**NO** Brukerveiledning

**SV** Bruksanvisning

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**SK** Návod na používanie

**HR** Upute za uporabu

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**HYDEAL II – 1P120(-KD) / 1P121(-KD)**  
*Instructions for use for Prosthetists/Orthotists*  
**Please read before using or working on the device**

1P12099  
2021-04

**Pass on § 3, 7, 8 and 9 of these instructions to the Patient**

**1. COMPONENTS INCLUDED**

Designation	Ref.	Included/Sold separately
Knee	1P120(-KD) or 1P121(-KD)	Included
One-part cosmetic cover	1G16	Sold separately
Short cylinder kit	1P11067	Sold separately

**2. DESCRIPTION, PROPERTIES AND MECHANISM OF ACTION**

**A. Description**


A 4-axis hydraulic knee available in four different versions:

- 1P120/1P121: Male pyramid proximal adapter, suitable for patients who have had a transfemoral amputation or hip disarticulation.
- 1P120-KD/1P121-KD: Threaded proximal adapter (M36×1.5), suitable for patients with a long amputation stump or who have had a knee disarticulation.

For all 4 versions, the distal connection is made via a Ø34 mm tube

**B. Properties**

Ref.	1P120	1P121	1P120-KD	1P121-KD
Colour	Grey	Orange	Grey	Orange
Weight	1,050 g		1,080 g	
Maximum flexion	145°			
Maximum patient weight (Including carried load)	125 kg			

**\*\* Dimensions in brackets are the dimensions when a short extension assist cylinder 1P11067 is used**

This device has undergone a three-million-cycle test, in accordance with standard NF EN ISO 10328 at a P6 load level (i.e. 125 kg), which corresponds to a service life of 4 to 5 years, depending on the patient's mobility level.

### C. Mechanism of action

When in extension, the system of linkages ensures that the instantaneous centre of rotation is to the rear of the load line, which keeps the knee stable during the stance phase. The speed of every flexion and extension movement can be adjusted independently.

An adjustment on the extension assist cylinder modifies the force applied to extend the knee.

Comment: a "Short cylinder" kit is available. It reduces the total height of the knee by 22 mm.

## 3. INTENDED USERS/INDICATIONS

This medical device is supplied to healthcare professionals (prosthetists/orthotists) who then instruct the patient in how to use it. The prescription is drawn up by a doctor. The doctor works with a prosthetist/orthotist to assess whether the patient is suited to using the device.

 This device should be used on ONE PATIENT ONLY. It must not be reused on other patients.

This device is exclusively intended to be used as prosthetic equipment for patients who have undergone transfemoral amputation, hip disarticulation or knee disarticulation. It is specifically recommended for active to very active patients (L3/L4) as it enables them to perform everyday activities as well as certain sports (cycling, running, rowing).

Maximum weight (including carried load): 125 kg

 **Not suitable for children.**

 The maximum knee flexion is 145°. This flexion may, however, be limited by the volume of the socket or by the cosmetic covering.

## 4. CLINICAL BENEFITS

The device allows:

- Everyday activities and certain sports (cycling, running, rowing) thanks to the absence of a cusp (or "over-centre" point) in flexion and a flexion angle of up to 145°.
- Changes in pace of walking, accommodated by the hydraulic unit.
- High levels of stability and safety in the stance phase thanks to the 4-axis joint geometry.

## 5. ACCESSORIES AND COMPATIBILITIES


Ref.		1P120 / 1P121	1P120-KD / 1P121-KD
Proximal adapter	Anchor	1K179 + 1K183 1K173 or 1K176 1K177 (max. weight 100 kg)	1K179
	Connectors	1K172 or 1K207-HD	1K184
Distal adapter	Tube and connector	1G01-HD or 1D52-P6	

## 6. ASSEMBLY AND FITTING TO THE PATIENT

### A. Assembly

 When tightening, comply with the following sequence and torques for the screws on the 34 mm diameter tube:

1. Tighten screw V1 to a torque of 5 Nm
2. Tighten screw V2 to a torque of 5 Nm
3. Retighten screw V1 to a torque of 5 Nm

#### Assembling the short cylinder kit (1P11067):

##### **Using the kit (sold separately):**

Suitable when a smaller overall knee size and/or less powerful extension assist is preferred.

##### **Installation:**

Remove the cylinder and the spring mounted on the knee and replace them with the elements supplied.

 Do not tighten the assembly forcefully until it locks.

##### **Guideline information:**

The new assembly at its maximum "screwed-in" setting applies a force equivalent to the original assembly tightened to 6 turns less than its maximum setting. The knee is then about 30 mm shorter overall.

### B. Aligning

#### Static alignment:

In the sagittal plane, the load line passes through the greater trochanter and between 0 and 10 mm in front of the axis of the knee.

Align the foot in compliance with the manufacturer's instructions.

In the frontal plane, the load line passes through the middle of the socket, through the middle of the knee and between the first and second toes.

**Dynamic alignment:**


Advancing the knee reduces the stability of the assembly and facilitates the initiation of the swing phase when walking.

Retracting the knee increases the stability and the feeling of safety for the patient when walking, but makes it more difficult to initiate flexion.

The positioning of the connector on the socket is important. If necessary, use a translating adapter to try out different position when determining the optimum alignment.

The threaded M36×1.5 adapter for the KD version does not accommodate alignment adjustment. The socket must therefore be made accordingly.

Check the alignment using suitable equipment (line laser, plumb line, etc.)

## C. Adjusting the walking pace

On delivery, the knee is set up for walking at medium pace (hex set screws E and F in their mid position, and extension assist cylinder R screwed in to its maximum minus three turns).

For patient safety reasons, the prosthetist/orthotist must begin the adjustment process with the knee set to its factory settings and with the patient walking between parallel bars.

- ⚠ The model of foot used may have an effect on the adjustment of the knee. The knee should be readjusted whenever the foot is changed. Prioritise patient safety in all circumstances, particularly to ensure that slopes and stairs can be negotiated safely.
- ⚠ After performing any maintenance the knee must be readjusted.

To adapt the knee to the patient's pace of walking, adjust hex set screws E and F using an Allen (hex) key.

It is advisable to perform these adjustments in the following order:

### 1. Adjusting the Flexion (F)

Turn hex set screw F: Screw in the set screw to reduce the angle of flexion, and unscrew to increase it.

*Adjustment range: 1 turn in the screw-in direction and 1 turn in the unscrew direction*

### 2. Adjusting the Extension (E)

Turn hex set screw E: Screw in to cushion the impact when the knee reaches full extension.


*Adjustment range: 1 turn in the screw-in direction and 1 turn in the unscrew direction*

### 3. Adjusting the extension assist (R)

If the optimal adjustment for the patient's walking pace cannot be achieved:

Adjust the extension assist cylinder R to change the assist power: Screw in for more powerful assist (faster extension) or unscrew for less powerful assist (slower extension).

- ⚠ After adjusting the assist power, the Extension (E) setting must be checked.

- ⚠ Once these adjustments have been made, ensure that the knee returns to full extension when walking slowly to ensure that there is no risk of the patient falling.

- ⚠ This hydraulic knee has been primarily designed to allow control during the swing phase. If, for example, the patient's body weight bears on the hydraulic system with the resistance to flexion set to high, there is a danger of degrading the knee.

## D. Finishing

It is advisable to reinforce the cosmetic cover in front of the joint to extend its service life.


### Maximum flexion stop:

When fully flexed, the socket (E) must never come into contact with the linkages (A) and (B).

The socket must be restrained by the maximum flexion stop (J) to preserve the knee.

- ⚠ The cosmetic cover may modify the initial knee settings and limit certain movements. It is advisable to always check the operation of the knee, once the cosmetic cover has been fitted.

- ⚠ Do not remove any of the covers over the knee.

## 7. DETECTING ANY MALFUNCTIONING

- ⚠ If you notice that the device is behaving abnormally or if its characteristics feel different, or if it has been badly knocked, consult your prosthetist/orthotist.

## 8. WARNINGS, CONTRAINDICATIONS, SIDE EFFECTS

### A. Warnings

- ⚠ Caution, the patient may need to take time to get used to a new knee.
- ⚠ To reduce the risk of falling when going down stairs, it is advisable to hold on to handrail.
- ⚠ So as not to damage the knee, do not apply talcum powder to eliminate rubbing noises - use silicon spray instead. Talc degrades the mechanical components, which might cause a malfunction and possibly a fall.  
**PROTEOR disclaims all liability if talcum powder is used.**
- ⚠ It is possible to trap a finger or pinch clothing in the knee joint. To avoid the any injury caused by the movement of the joint, make sure that no-one's hand or fingers is near the mechanism.
- ⚠ **The knee is resistant to bad weather, but needs to be dried if it does get wet.**

### B. Contraindications

- ⚠ The following are strictly prohibited: removing the covers, screwing in or unscrewing any screws on the knee with the exception of the two adjustment screws (**E** and **F**) and the two tightening screws (**V1** and **V2**) on the Ø34 tube clamp adapter.
- ⚠ It must be possible to rotate the extension assist cylinder R. Never tighten it fully so that it can be loosened if necessary.
- ⚠ Never lubricate the knee pins; this can cause them to degrade rapidly.
- ⚠ The knee has been designed for patients with a maximum weight of 125 kg (including carried load). If overloaded, a safety feature releases the hydraulic system which can lead to rapid flexion of the knee.
- ⚠ If a load is being carried, the operation of the knee may be disrupted. Particularly, the resistance when going down stairs may be insufficient, or the knee may lock up suddenly.
- ⚠ If the patient's weight increases significantly, he or she should arrange for the prosthetist/orthotist to readjust the knee's safety settings.
- ⚠ The warranty does not cover damage or degradation caused by misuse, unsuitable alignment, use without suitable protection in a very dusty environment or any other inappropriate use.
- ⚠ The knee should not be exposed to environments which might induce corrosion of the metal components (fresh water, sea water, chlorinated water, acids, etc.).
- ⚠ **Showering or taking a bath while wearing the prosthesis is prohibited;** this might impair its properties and its correct operation.
- ⚠ Intensive use of the knee brake (e.g. due to a prolonged period of walking down a slope or stairs) may cause the hydraulic system to overheat and the braking effect to diminish. Do not touch the knee and reduce the activity level to give the components time to cool down.
- ⚠ Never leave this device near a source of heat since it might burn or release toxic fumes.
- ⚠ If used in very cold (<10°C) or very hot (>40°C) conditions the behaviour of the knee may change significantly. In this case, take precautions when walking and going down slopes or stairs.
- ⚠ The use of solvents is prohibited.

### C. Side effects

There are no side effects directly associated with the device.

Any serious incident that occurs which relates to the device must be reported to the manufacturer and to the competent authority of the member state in which the incident occurred.

## 9. MAINTENANCE, STORAGE, DISPOSAL AND SERVICE LIFE

### A. Maintenance/cleaning

- ⚠ You can clean the knee using a damp sponge
- ⚠ Do not immerse the knee or hold it under running water
- ⚠ Dry you knee if it is exposed to bad weather (rain) or accidental splashing.

The following components may need to be replaced:

- Spare extension stop: EE015 (x1) to be bonded using PROTEOR XC051 glue
- Spare maximum flexion stop: EY001 (x1), self-adhesive

### B. Storage

- ⚠ Service and storage temperature: -10°C to +40°C

Relative humidity of the air: no restrictions

### C. Disposal

The various components of this device are considered as special waste: elastomer, plastic, aluminium, titanium, steel, brass and oil. They must be treated in accordance with current, applicable legislation.

### D. Service life

It is advisable to arrange for a prosthetist/orthotist to check the product once a year.

## 10. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of first declaration of conformity
---	--------------	---	-----------------	---	--

## 11. REGULATORY INFORMATION

This product is a CE-marked medical device that is certified as conforming with regulation (EU) 2017/745

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PROTEOR SAS