HYTREK

1P131 / 1P131-KD

PROFESSIONAL

FR	Notice d'Utilisation	FI	Käyttöohjeet
EN	Instruction for use	PL	Instrukcja użytkowania
DE	Gebrauchsanweisung	CS	Návod k použití
IT	Istruzioni per l'uso	SK	Návod na používanie
ES	Instrucciones de uso	HR	Upute za uporabu
PT	Instruções de utilização	RU	Инструкция по использованию
NL	Gebruikershandleiding	UK	Інструкція з використання
DA	Brugervejledning	JA	取扱説明書
NO	Brukerveiledning	ZH	使用说明
SV	Bruksanvisning	AR	تعليمات الاستخدام

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HYTREK – 1P131(-KD)
Instructions for use for Prosthetists/Orthotists
Please read before using or working on the device

1P13099
2022-11

It is important to read and explain the patient information leaflet to the patient.

1. COMPONENTS INCLUDED

Designation	Ref.	Included/Sold separately
Knee	1P131(-KD)	Included
Protective hood	1P13055	Included
2.5 mm Allen (hex) key	WKF026	Included
Cover	1P13029-N	Included
	1P13028-N	
	1P13027-N	
One-part cosmetic cover	1G16	Sold separately

2. DESCRIPTION, PROPERTIES AND MECHANISM OF ACTION

A. Description

A single-axis knee offering hydraulic regulation of the stance and swing phases, fitted with a manual lock that halts the flexion of the knee.


It is available in four different versions:

- 1P131: Male pyramid proximal connector, suitable for patients who have had a transfemoral amputation or hip disarticulation
- 1P131-KD: Threaded proximal connector (M36×1.5), suitable for patients with a long transfemoral amputation stump or who have had a knee disarticulation

For both versions, the distal connection is made via a $\varnothing 34$ mm tube

B. Properties

Ref.	1P131	1P131-KD
Weight	1265 g	1285 g
Maximum flexion	120°	
Maximum patient weight (Including carried load)	150 kg	

This device has undergone a three-million-cycle test, in accordance with standard NF EN ISO 10328 at a P7 load level (i.e. 150 kg), which corresponds to a service life of 4 to 5 years, depending on the patient's mobility level.

C. Mechanism of action

The knee contains an oil-filled hydraulic cylinder.

The hydraulic system and the control of its operation provide safe management of the stance phase (transition after the swing phase)

and resistance to flexion) and the swing phase (impact at the end of extension and dampening of flexion). The angle of flexion of the knee is limited mechanically to 120°. This amplitude lets the patient kneel down and ride a bicycle.

It offers four separate adjustments, extension assist and a system for locking the flexion.

The knee’s intrinsic characteristics (flexion, braking, hydraulic regulation) allow patients to negotiate slopes and stairs using alternating feet. As a precautionary measure, it is advisable to hold on to the handrail when walking down stairs.

The braking torque when descending stairs is high (150 Nm) up to a significant flexion angle (75°). This torque gradually reduces as the flexion angle reduces.


Simply by pressing on a flexion locking button, the knee can be locked to prevent any flexion. With regard to extension, the lock halts the flexion and allows the knee to extend. The patient can thus lock the knee to make it easier to negotiate obstacles, stand without tiring and perform activities safely (e.g. climbing up a ladder).

Pressing a second time on the button unlocks the knee.

3. INTENDED USERS/INDICATIONS

This medical device is supplied to healthcare professionals (prosthetists/orthotists) who then instruct the patient in how to use it. The prescription is drawn up by a doctor. The doctor works with a prosthetist/orthotist to assess whether the patient is suited to using the device.

⚠ This device should be used on ONE PATIENT ONLY. It must not be reused on other patients.

This device is exclusively intended to be used as prosthetic equipment for patients who have undergone transfemoral amputation, hip disarticulation or knee disarticulation. It is specifically recommended for active to very active patients (L3/L4) as it allows them to carry out their everyday activities on all types of terrain as well as negotiate slopes and stairs.

Maximum weight (including carried load): 150 kg

⚠ **Not suitable for children.**

⚠ The maximum flexion of the knee is 120°. This flexion may, however, be limited by the volume of the socket or by the cosmetic covering.

4. CLINICAL BENEFITS

The device allows:

- The patient to walk down stairs and slopes by bearing weight on the amputated leg to limit the effort on the sound leg.
- The knee to be locked in flexion using a manual lock to allow certain activities to be performed safely (e.g. climbing a ladder).
- The prosthetist/orthotist to adjust certain settings to adapt the knee to the patient.
- Walking at different speeds thanks to the hydraulic regulation.
- The hip to do less work since some of the knee extension force is provided by an extension assist spring.

5. ACCESSORIES AND COMPATIBILITIES


Ref.		1P131	1P131-KD
Proximal adapter	Anchor	1K179+ 1K183 (max. weight 150 kg) 1K173 or 1K176 (max. weight 125 kg) 1K177 (max. weight 100 kg)	1K179 (max. weight 150 kg)
	Connectors	1K172 or 1K207-HD (max. weight 150 kg)	1K184 (max. weight 150 kg)
Distal adapter	Tube and connector	1G01-HD (max. weight 150 kg) 1D52-P6 (max. weight 125 kg)	

6. ASSEMBLY AND FITTING TO THE PATIENT

A. Assembly

⚠ When tightening, comply with the following sequence and torque for the screws on the 34 mm diameter tube:

1. Tighten screw V1 to a torque of 5 Nm
2. Tighten screw V2 to a torque of 5 Nm
3. Retighten screw V1 to a torque of 5 Nm

B. Aligning

Static alignment:

In the sagittal plane, the load line passes through the greater trochanter and between 0 and 5 mm in front of the axis of the knee.

Ensure that any flexion contracture in the patient is accommodated.

Align the foot in compliance with the manufacturer's instructions.

In the frontal plane, the load line passes through the middle of the socket, through the middle of the knee and between the first and second toes.

Dynamic alignment:


Advancing the knee reduces the stability of the assembly and facilitates the initiation of the swing phase when walking.

Retracting the knee increases the stability and the feeling of safety for the patient when walking, but makes it more difficult to initiate flexion.

The positioning of the connector on the socket is important. If necessary, use a translating adapter to try out different position when determining the optimum alignment.

The threaded M36x1.5 adapter for the KD version does not accommodate any alignment adjustment. The socket must therefore be made accordingly.

Check the alignment using suitable equipment (line laser, plumb line, etc.)


C. Adjustment

For patient safety reasons, the prosthetist/orthotist must begin the adjustment process with the knee set to its factory settings and with the patient walking between parallel bars.

- ⚠ The model of foot used may have an impact on the adjustment of the knee. The knee should be readjusted whenever the foot is changed.
- ⚠ Prioritise patient safety in all circumstances, particularly to ensure that slopes and stairs can be negotiated safely. After performing any maintenance the knee must be readjusted.

A 2.5 mm Allen (hex) key, supplied with the knee, is used to make all the adjustments.

Adjust as described below, ideally in **the order indicated**:

1. Adjusting the intensity of the resistance in the stance phase (Screw R)

- ⚠ The HYTREK knee is delivered with a fairly high resistance setting. The initial adjustments should be performed with care and with safety in mind to avoid any risk of the patient falling. Screw in the screw, 1/12th of a turn at a time, to increase the resistance in the stance phase (braking). Unscrew to produce the opposite effect. Then adjust more finely if required.
Adjustment range = ¾ of a turn (1 turn = same adjustment)
With the patient between parallel bars, and in order to define an initial stance resistance adjustment, ask the patient to sit down on a chair. Then, as a follow-up step after performing all the other adjustments on the knee, refine the resistance with the patient walking, with alternating feet, down stairs and then down a slope, and always with a handrail next to the patient for safety reasons.

2. Adjusting the threshold value for the transition to swing phase (Screw S)

- ⚠ The HYTREK knee is delivered with this adjustment set to a value that is suitable for most users. The dynamic alignment should allow the transition to swing phase and the activation of braking. The threshold may be modified if the swing phase does not initiate at the instant when the toes leave the ground (toe-off), or if resistance is not in place at heel strike. Screw in, one ¼ of a turn at a time, until there is no locking of the knee at toe-off. Unscrew, one ¼ of a turn at a time, to facilitate the triggering of the brake, then fine adjust one 1/12th of a turn at a time.

Adjustment range = -1 turn/+2 turns

- ⚠ Do not fully unscrew screw S since this might damage the mechanism. If in doubt, reset the adjustment by turning the screw until its head is flush with the surface.
- ⚠ There is a significant risk of the patient falling if this adjustment is set incorrectly.


3. Adjusting the flexion during the swing phase (Screw F)

Screw in, one 1/12th of a turn at a time, to limit the flexion during the swing phase, notably at fast walking paces. Unscrew to produce the opposite effect. Adjust more finely if required.
Adjustment range: ¾ of a turn (1 turn = same adjustment)

4. Adjusting the absorption of the terminal impact (Screw E)

Screw in, one ¼ of a turn at a time, to increase the absorption of the terminal impact (at the end of extension). Unscrew to produce the opposite effect, and then adjust more finely if required.
Adjustment range: ± 1 turn

- ⚠ Once these adjustments have been made, ensure that the knee returns to full extension when walking slowly to ensure that there is no risk of the patient falling.

Maximum flexion stop

When fully flexed, the socket may come into contact with the posterior band, located above the adjustment buttons.

The socket must not come into contact with the hydraulic system since this might damage or degrade the knee.

The knee is fitted with an internal stop which limits its flexion to 120°.

Lock

The HYTREK knee is fitted with a lock which can be activated when the patient is standing or sitting.

To lock the knee, press on the button at the rear of the knee. Always check that the lock has been engaged before applying load to the knee. Press on the button next to the lock button to unlock the knee.

D. Finishing

It is advisable to reinforce the cosmetic cover in front of the joint to extend its service life.


If worn without a cosmetic cover, the protective hood supplied may be affixed to the knee, using adhesive (XC050).

The cosmetic cover may modify the knee's initial settings and limit certain movements. It is advisable to always check the operation of the knee once the cosmetic cover has been fitted.

Do not remove any of the knee's covers, notably the posterior band which protects the hydraulic system when the knee is fully flexed.

7. DETECTING ANY MALFUNCTIONING

If you notice that the device is behaving abnormally or if its characteristics feel different, or if it has been badly knocked, consult your prosthetist/orthotist.

8. WARNINGS, CONTRAINDICATIONS, SIDE EFFECTS

A. Warnings

Caution, a knee with an electronic braking mechanism does not behave like a braked knee; the patient may need some time to get used to a new knee when switching between the two.

To reduce the risk of falling when going down stairs, it is advisable to hold on to handrail.

If the knee is not used for a few hours, it may make a clicking or cracking noise when first used. This does not impair its operation.

So as not to damage the knee, do not apply talcum powder to eliminate rubbing noises - use silicon spray instead. Talc degrades the mechanical components, which might cause a malfunction and possibly a fall.

PROTEOR disclaims all liability if talcum powder is used.

It is possible to trap a finger or pinch clothing in the knee joint. To avoid the any injury caused by the movement of the joint, make sure that no-one's hand or fingers is near the mechanism.

If the knee is immobilised for a prolonged period in very cold conditions (<10°C), the hydraulic system will not return to normal operation until the patient has completed a few steps. Walk carefully until it feels normal again.

The knee is resistant to bad weather, but needs to be dried if it does get wet.

B. Contraindications

The following are strictly prohibited: removing the covers, screwing in or unscrewing any screws on the knee with the exception of the four adjustment screws (**R**, **S**, **F** and **E**) and the two tightening screws (**V1** and **V2**) on the Ø34 tube clamp adapter.

Never lubricate the knee pins; this might cause them to degrade rapidly.

The knee has been designed for patients with a maximum weight of 150 kg (including carried load). If overloaded, a safety feature releases the hydraulic system which can lead to rapid flexion of the knee.

If a load is being carried, the operation of the knee may be disrupted. Particularly, the resistance when going down stairs may be insufficient, or the knee may lock up suddenly.

If the patient's weight increases significantly, he or she should arrange for the prosthetist/orthotist to readjust the knee's safety settings.

- ⚠ The warranty does not cover damage or degradation caused by misuse, unsuitable alignment, use without suitable protection in a very dusty environment or any other inappropriate use.
- ⚠ The knee should not be exposed to environments which might induce corrosion of the metal components (fresh water, sea water, chlorinated water, acids, etc.).
- ⚠ **Showering or taking a bath while wearing the prosthesis is prohibited;** this might impair its properties and its correct operation.
- ⚠ Intensive use of the knee brake (e.g. due to a prolonged period of walking down a slope or stairs) may cause the hydraulic system to overheat and the braking effect to diminish. Do not touch the knee and reduce the activity level to give the components time to cool down.
- ⚠ Never leave this device near a source of heat since it might burn or release toxic fumes.
- ⚠ If used in very cold (<10°C) or very hot (>40°C) conditions the behaviour of the knee may change significantly. In this case, take precautions when walking and going down slopes or stairs.
- ⚠ The use of solvents is prohibited.

C. Side effects

There are no side effects directly associated with the device.

Any serious incident that occurs which relates to the device must be reported to the manufacturer and to the competent authority of the member state in which the incident occurred.

9. MAINTENANCE, STORAGE, DISPOSAL AND SERVICE LIFE

A. Maintenance/cleaning

- ⚠ You can clean the knee using a damp sponge
- ⚠ Do not immerse the knee or hold it under running water
- ⚠ Dry you knee if it is exposed to bad weather (rain) or accidental splashing.

B. Storage

- ⚠ Service and storage temperature: -10°C to +40°C
- Relative humidity of the air: no restrictions

C. Disposal

The various components of this device are considered as special waste: elastomer, plastic, aluminium, titanium, steel, brass and oil. They must be treated in accordance with current, applicable legislation.

D. Service life

It is advisable to arrange for a prosthetist/orthotist to check the product once a year.

10. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of first declaration of conformity
---	--------------	---	-----------------	---	--

11. REGULATORY INFORMATION

This product is a CE-marked medical device that is certified as conforming with regulation (EU) 2017/745

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PROTEOR SAS