ShockWave™

F20

- | | | | |
|-----------|--------------------------|-----------|-----------------------------|
| FR | Notice d'Utilisation | FI | Käyttöohjeet |
| EN | Instruction for use | PL | Instrukcja użytkowania |
| DE | Gebrauchsanweisung | CS | Návod k použití |
| IT | Istruzioni per l'uso | SK | Návod na používanie |
| ES | Instrucciones de uso | HR | Upute za uporabu |
| PT | Instruções de utilização | RU | Инструкция по использованию |
| NL | Gebruikershandleiding | UK | Інструкція з використання |
| DA | Brugervejledning | JA | 取扱説明書 |
| NO | Brukerveiledning | ZH | 使用说明 |
| SV | Bruksanvisning | AR | تعليمات الاستخدام |

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ShockWave™ – F20

Instructions for use for prosthetists

Read before use

IFU-01-097
Rev. F
2023-10

Pass on § 12 to 18 of these instructions to the patient.

1. INCLUDED ITEMS

Part description	Part number	Included / Sold separately
ShockWave Foot	F20-S3-xxAxx-xx*	Included
Foot Shell	FTC-3F-xxxxx-Rx* FTC-2F-xxxxx-Sx*	Sold separately
Black spectra sock	S0-NPS-200xx-00*	Suitable sock included
Stiffening bumpers	KIT-00-1146U-00	Included


* see catalog

2. DESCRIPTION

ShockWave™ is an energy-return prosthetic foot consisting of:

- A J-shaped carbon keel
- A carbon sole plate
- A vertical loading pylon that provides controlled motion for shock absorption and axial rotation
- A spectra sock

Available with regular toe or right or left sandal toe and delivered with a set of elastomer stiffening bumpers.

3. PROPERTIES

Side	Right or Left		
Size	23 to 25 cm	26 to 28 cm	29 to 30 cm
Sandal toe option	Yes	Yes	No
Weight*	897 g / 1.98 lb	978 g / 2.16 lbs	1131 g / 2.49 lbs
Build height*	176 mm / 6.9"	189 mm / 7.4"	200 mm / 7.9"
Heel height	10 mm / 3/8"		

*Based on a sizes 23, 26, 29, cat. 4, with Foot Shell, Spectra Sock and 10 mm heel height

This device has been tested according to ISO 10328 for a maximum patient weight up to 147 kg for 2 million cycles.

Selection of foot category based on patient's weight and activity level										
Weight ^{*)}	lb	110-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325
	kg	44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
Activity level	Low	1	1	2	3	4	5	6	7	8
	Moderate	1	2	3	4	5	6	7	8	-
	High	2	3	4	5	6	7	8	-	-

^{*)} Body mass limit not to be exceeded (ISO 10328)

4. MECHANISM OF ACTION


From Initial Contact and into Loading Response, the sole and heel store energy, which is released as the user rolls into midstance. As the user surpasses midstance the forefoot is loaded, storing energy that is released at terminal stance for propulsion forward. The shock unit cushions shock and dampens rotational forces throughout the gait cycle and then returns the foot to a neutral position.

5. INTENDED USE/INDICATIONS

This medical device is supplied to healthcare professionals (prosthetists) who will train the patient in its use. The prescription is drawn up by a doctor who assess the patient's ability to use the device.

⚠ This device is for multiple use on a **SINGLE PATIENT**. It must not be used on another patient.

This device is intended to be integrated in a custom-made external lower limb prosthesis to ensure the function of the foot in patients with unilateral or bilateral lower limb amputation and/or congenital limb deficiencies.

This device is indicated for patients with moderate to high activity level (K3 to K4) for walking and physical activity without excessive overload.

Maximum weight (load carrying included): 147 kg / 325 lb (See table §3)

6. CLINICAL BENEFITS

- Walking comfort
- Possible walking on uneven ground
- Reduction of socket forces

7. ACCESSORIES AND COMPATIBILITY

An appropriate foot shell must be mounted on the foot module (refer to our catalog).


The foot includes a male pyramid connection designed to be compatible with standard female pyramid connectors (refer to our catalog).

8. ALIGNMENTS

Bench Alignment

Before installing the prosthesis on the patient:

- Align the foot module in plantar/dorsal flexion, by inserting the foot (with foot shell) into the shoe to account for the heel height
- Align the socket in adduction/abduction, to ensure an appropriate angle in the frontal plane
- Align the socket in flexion/extension, to ensure an appropriate angle in the sagittal plane
- Adjust the position of the socket so that the load line falls through the center of the pylon (see illustration)

Dynamic Alignment

To optimize roll-over from heel to forefoot, adjust the following variables:

- Foot position in the anterior/posterior plane
- Plantar/dorsal flexion
- Heel flexibility

The dynamic alignment is performed in accordance with good professional practices.

9. ASSEMBLY

The ShockWave is pre-assembled and consists of a foot module, a Spectra sock and a foot shell. After dynamic alignment, tighten the pyramid adjustment screws according to the specifications of the connector manufacturer. Secure pyramid adjustment screws with a thread locking adhesive (i.e., Loctite 242).

Spectra Sock

A Spectra sock is included to protect the foot shell and minimize noise. It must be placed over the foot module before mounting the foot shell.

Foot Shell

To install and remove the foot shell, use the Foot Shell Removal Tool to prevent damage to the foot module.

- ⚠ Never remove the foot from the foot shell by pulling manually. Never use a screwdriver or any other inappropriate instrument to remove it. This could damage the foot.

10. ADJUSTMENTS

The heel stiffness at foot loading can be adjusted using the stiffening bumpers. The bumpers may be temporarily attached using the pre-applied adhesive in the location indicated in the table below. The recommended location will stiffen the heel of about 1

category. If necessary, move the bumpers anteriorly (stiffer) or posteriorly (softer) to get the desired stiffness. For permanent placement, clean off the pre-applied adhesive with acetone and adhere bumpers using cyanoacrylate glue.

Foot size (A)	Distance from the rear end of the sole (B)		
23-25 cm	29 mm / 1.1"		
26-28 cm	35 mm / 1.4"		
29-30 cm	48 mm / 1.9"		

11. TROUBLESHOOTING

CONCERN	SYMPTOM	SOLUTION
Heel too soft	<ul style="list-style-type: none"> • Foot flat occurs too rapidly • Toe feels excessively stiff • Knee hyperextension 	<ul style="list-style-type: none"> • Shift socket anteriorly in relation to the foot • Attach stiffening bumpers. See section 10 above for installation details
Heel too hard	<ul style="list-style-type: none"> • Rapid knee flexion, instability • Heel to toe progression too fast • Lack of energy return sensation 	<ul style="list-style-type: none"> • Shift socket posteriorly in relation to the foot • Verify appropriate foot module category
Foot module too stiff	<ul style="list-style-type: none"> • Flat spot in rollover motion at slow cadences 	<ul style="list-style-type: none"> • Consider a lower category foot module
Foot module too soft	<ul style="list-style-type: none"> • Clicking noise at initial contact • Excessive toe deflection with high impact activity 	<ul style="list-style-type: none"> • Consider a higher category foot module

12. WARNINGS

- ⚠ In case of damaged packaging, check the integrity of the device.
- ⚠ Never use the foot without a foot shell and a Spectra sock.
- ⚠ Never loosen the pyramid fastening screws.
- ⚠ The patient must inform their prosthetist if they gain or lose significant weight.
- ⚠ Always use the foot with a sock and shoe. Failure to adhere to this advice may cause product failure, as well as serious injury.
- ⚠ Make sure that the foot and inside of the foot shell are free of impurities (e.g., sand). The presence of impurities causes the graphite parts to wear out. Clean the foot according to the instructions (see §16).
- ⚠ After swimming, using in water or if splashed by a liquid, the foot must be cleaned (see §16).
- ⚠ If the patient notices any abnormal behavior or feel any changes in the characteristics of the device (noise, play, excessive wear...), or if the device has received a severe impact, they should stop using the device and consult their prosthetist.
- ⚠ Failure to follow the instructions for use is dangerous and will void the warranty.

13. CONTRAINDICATIONS

- ⚠ Do not use for a patient whose maximum weight (load carrying included) may exceed 147 kg / 325 lbs.
- ⚠ Do not use for activities associated with a risk of significant impact or excessive overloading.
- ⚠ Do not use for a patient who do not meet the requirements of a Medicare functional level of K3 or higher.

14. SIDE EFFECTS

There are no known side effects directly associated with the device.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user is established.

15. MAINTENANCE AND CONTROL

The foot module must be inspected by the prosthetist at least every six months. Inspections at shorter intervals are required if the user is more active.

The Spectra sock and the foot shell must be replaced by the prosthetist at regular intervals, depending on the patient's level of activity. If these parts are damaged, it can lead to premature foot wear.

The lifetime of the foot depends upon the patient's level of activity.

16. CLEANING

Remove the foot shell and the Spectra sock, clean the foot with clear water and neutral soap and dry carefully.

The foot shell can be cleaned with a damp cloth or sponge.

 The device is not resistant to solvents. Exposure to solvents may cause damage.

17. ENVIRONMENTAL CONDITIONS

Temperature range for use and storage: -29 to 49°C [-20 to 120° F]

Relative air humidity: no restrictions

Water-resistant: the device is resistant to occasional immersion and splashes.

 The foot shell is not resistant to ultraviolet light (UV). Do not store it under direct sunlight.

18. DISPOSAL

The foot is made of carbon fiber composite, epoxy, titanium, stainless steel, and rubber. The device and its packaging must be disposed of in accordance with local or national environmental regulations.

19. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of 1 st declaration
	Authorized representative in the European Union		Single patient, multiple use		

20. REGULATORY INFORMATION

This product is a CE-marked medical device and is certified as conforming with Regulation (EU) 2017/745.

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PROTEOR SAS