EXONEO CORE®

05CxCx-Px

EN	Instruction for use	FI	Käyttöohjeet
FR	Notice d'Utilisation	PL	Instrukcja użytkowania
DE	Gebrauchsanweisung	CS	Návod k použití
IT	Istruzioni per l'uso	SK	Návod na používanie
ES	Instrucciones de uso	HR	Upute za uporabu
PT	Instruções de utilização	RU	Инструкция по использованию
NL	Gebruikershandleiding	UK	Інструкція з використання
DA	Brugervejledning	JA	取扱説明書
NO	Brukerveiledning	ZH	使用说明
SV	Bruksanvisning	AR	تعليمات الاستخدام

	EXONEO CORE® <i>Instructions for use</i> Please read before use	IFU-02-022 Rev. 1 2026-02
---	---	---

Provide the patient with points 10, 11, 12, 13, 14, 15, 16 and 17 of these instructions.

1. COMPONENTS INCLUDED

Name	Reference	Included / Sold separately
EXONEO CORE prosthesis	05CxCx-Px-xx*	Included
Foot shell assembly	Large (size 25 to 29): 0100x01* Small (size 21 to 24): 0400x01*	Included

*see catalogue

2. DESCRIPTION

EXONEO CORE is an energy-return prosthetic foot composed of:

- A fibreglass blade
- One or more tension springs
- A shock-absorbing component in the heel
- A forefoot
- A male pyramid connector
- A protective foot shell

EXONEO CORE is delivered fully assembled. The foot shell must be fitted onto the foot by a prosthetist.

3. PROPERTIES

Type of prosthesis	Small size	Large size
Size	21 to 24 cm	25 to 29 cm
Side	Right or Left	Right or Left
Weight	403 to 465 g	593 to 739 g
Patient weight	40 to 90 kg	51 to 150 kg
Build height	65 mm	75 mm
Total height	80 mm	90 mm

This device has been tested as per standard ISO 10328 for a loading level of P2 to P7 for 2 million cycles, corresponding to a service life of 3 years, depending on the patient's activity level.

Selection of foot category according to patient weight and foot size:							
Weight ^{*)}	in kg	40-50	51-70	71-90	91-110	111-130	131-150
Foot size	21-24 cm	X	X	X			
	25-29 cm		X	X	X	X	X
Foot configuration (and loading level tested)		P2	P3	P4	P5	P6	P7

^{*)} Body mass should not exceed this limit (ISO 10328)

4. ACTION MECHANISM

The main feature of this prosthetic foot is that it is polyarticulated to reproduce the anatomical movement of the human foot. As such, it is made up of three interconnected parts: the forefoot, the ankle and the heel.


During stance, the heel cushions the initial impact, then the blade deforms and the spring or springs are stretched to store energy, which is released at the end of stance.

5. INTENDED USERS/INDICATIONS

This medical device is supplied to healthcare professionals (prosthetists), who will train the patient in its use. The prescription is drawn up by a doctor, who assesses the patient's ability to use the device.

⚠ This device is for multiple use on a **SINGLE PATIENT**. It must not be used on another patient.

This device is designed to be incorporated into a tailor-made external lower-limb prosthesis to provide foot function for patients who are bilateral or unilateral amputees and/or have a congenital lower-limb defect.

This device is indicated for patients with a low activity level (CIF d4601) for non-sport walking (outdoors or indoors) and for going up and down stairs.

Maximum weight (including carried load): 150 kg / 330 lbs (See section 3)

Minimum weight: 40 kg / 88 lbs (See section 3)

Foot length: 21 to 29 cm (See section 3)

6. CLINICAL BENEFITS

The device provides everyday stability and comfort for patients with transfemoral or transtibial amputations walking at normal speeds over smooth or uneven ground.

7. ACCESSORIES AND COMPATIBILITY


The foot includes a male pyramid connector, designed to be compatible with standard female pyramid adapters (see our catalogue).

8. ALIGNMENTS

Static alignment

Align the device such that the load line passes through the posterior of the foot, 1/3rd of the length from the back.

Before fitting the prosthesis to the patient:

- EXONEO CORE has been designed without the usual 10 mm heel lift. Alignment can be done with or without shoes, but it is preferable to make adjustments with shoes on. Align the foot in plantar flexion/dorsiflexion by placing the foot with the foot shell in the shoe (if necessary) to take account of the height of the heel.
- Align the socket in adduction/abduction to ensure an appropriate angle in the frontal plane.
- Align the socket in flexion/extension to ensure an appropriate angle in the sagittal plane.

Dynamic alignment

To optimise foot roll-over, adjust the following variables:


- Anteroposterior translation
- Plantar flexion/dorsiflexion

Dynamic alignment is performed in accordance with good professional practices.

9. ASSEMBLY

⚠ **NB:** do not mix parts from different configurations (P2, P3, P4, P5, P6, P7).

ASSEMBLING THE PROTECTIVE FOOT SHELL (use a shoehorn):


REMOVING THE PROTECTIVE FOOT SHELL (use a shoehorn):

Insert the shoehorn at the back of the heel, between the foot shell and the prosthesis.

NB:

The shoehorn should be pushed slightly towards the foot shell to avoid breaking the bottom of the heel during insertion.

Once the heel is clear, the foot shell can be safely removed.

10. DETECTING MALFUNCTIONS

- ⚠ If the patient notices any abnormal behaviour or feels any changes in the characteristics of the device (noise, play, excessive wear, etc.), or if the device has received a severe impact, then it should no longer be used and the patient should consult their prosthetist.

11. WARNINGS

- ⚠ Keep the label on the box for the life of the prosthetic foot. In the event of a complaint, provide the information contained on the label.
- ⚠ Do not reuse the protective foot shell from another prosthetic foot, whether or not it is the same size.
- ⚠ Check that the device is complete if the packaging is damaged.
- ⚠ Bilateral patients – allow for an extended evaluation period: the initial feeling of a stable and comfortable heel may change over time.

Warning to be conveyed to the patient:

- ⚠ Do not jump, run or climb with this foot.
- ⚠ Never use the foot without the protective foot shell.
- ⚠ Do not remove the foot shell without a shoehorn.
- ⚠ The prosthetic foot and its foot shell must not be exposed to chemical substances (oil, fuels, acidic or basic solutions, etc.).
- ⚠ Do not leave in direct sunlight or behind glass. Do not store near a heat source (fireplace, radiator).
- ⚠ Check the integrity of the foot shell and/or components before use.
- ⚠ If you change your shoes and this significantly increases or decreases the height of the heel, please consult your prosthetist.
- ⚠ In the event of a significant change in weight, please consult your prosthetist.
- ⚠ Failure to respect these instructions for use is dangerous and will invalidate the warranty.

12. CONTRAINDICATIONS

- ⚠ Do not use for patients whose maximum weight (including carried load) exceeds 150 kg / 330 lbs.
- ⚠ Do not use for activities where there is a risk of significant impact or excessive overloading.

13. SIDE EFFECTS

This device has no known direct side effects.

Any serious incident that occurs relating to the device must be reported to the manufacturer, and the competent authority within the Member State where the user resides.

14. MAINTENANCE AND CHECKS

No maintenance operations such as lubrication, working on screws or other parts are required.

The foot must be inspected by a prosthetist at least once a year. Inspections at shorter intervals are required if the user is more active.

The protective foot shell must be replaced by the prosthetist at regular intervals, depending on the patient's level of activity.

If you experience any difficulties with your prostheses, please contact your prosthetist.

15. CLEANING THE DEVICE

The prosthetic foot and its foot shell are resistant to chlorinated, salty and fresh water. If exposed to these, immediate maintenance is required. Empty the protective foot shell when it contains water and rinse the whole unit with clean water.

The prosthetic foot and its protective foot shell must be cleaned regularly (at least once a week), and visually inspected for any abnormal signs of wear. For maintenance, use soapy water only, and wash by hand. Dry with a clean cloth or in the open air. Do not put in the oven, microwave or washing machine.

Traces of wear and tear (oxidation spots, discolouration, etc.) may appear with normal use and maintenance. These marks are purely aesthetic and do not alter the technical characteristics of the foot. In the event of pronounced signs of abnormal wear, a diagnosis by a prosthetist is necessary.

16. ENVIRONMENTAL CONDITIONS

Temperature range for use and storage: -20°C to +60°C

Relative air humidity: no restrictions

Water resistance: fresh, salt and chlorinated.

17. DISPOSAL

This device is made of: carbon fibre, stainless steel, galvanised steel, titanium, anodised aluminium and plastic. The device and its packaging must be disposed of in accordance with applicable local or national environmental regulations.

18. DESCRIPTION OF SYMBOLS

	Manufacturer		Identified risk		CE marking and year of first declaration		Single patient, multiple use
---	--------------	---	-----------------	---	--	---	------------------------------

19. REGULATORY INFORMATION

This device is a CE-marked medical device that is certified as conforming with Regulation (EU) 2017/745.

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PROTEOR SAS