RUSH FOOT®

by  PROTEOR

RUSH Foot®

Instructions for Use

07/2021

Table of Contents / Tabla de contenido / Sommaire / Inhaltsangabe Indice / Obsah / 目次

**General Product Statement / Información del producto /
Information produit / Produktinformationen /
Informazioni sul prodotto / Informace o výrobku / 製品情報**

Page 3

English
Pages 4 - 8

Spanish / Español
Pages 9 - 13

French / Français
Pages 14 - 18

German / Deutsch
Pages 19 - 23

Italian / Italiano
Pages 24 - 28

Czech / Český
Pages 29 - 33

Japanese / 日本語
Pages 34 - 38

General Product Statement

ENGLISH : The RUSH Foot[®]'s superior engineering delivers the most realistic and responsive foot and ankle motion available. It's unique glass composite material eliminates prevalent weak points and is exceedingly more flexible than carbon fiber. The RUSH Foot[®] delivers an unmatched range of motion from the roll-through design which allows for a smooth transition from heel strike through toe-off with no "dead spots" in between. The innovative design helps maintain a natural, smooth gait, with superior energy return, even in rugged or uneven terrain. The radically simplistic design demonstrates that stress points won't fail when subjected to pressure, water, or impact. Virtually indestructible even in the most extreme conditions.

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RUSH HiPro, RAMPAGE, & RAMPAGE LP

Instructions for Use

**This instructional booklet is intended for the Prosthetist and the patient.
Read it carefully and completely before any use.**

**Please transmit to the patient the safety instructions and all the information
necessary for the proper use of the product (Sections 8 & 9).**

1—Description and Functions

The fiberglass composite material of the RUSH HiPro and Rampage feet provide great flexibility, resistance, and durability.

Rampage LP feet offer the same great flexibility, resistance, and durability as the HiPro and RAMPAGE, in a lower build height configuration.

Each RUSH foot[®] is comprised of the following components:

- An upper blade (Instep) that provides energy storage and return at toe-off.
- A continuous lower blade (Half Rocker Sole Plate) with a patented design that ensures progressive flattening of the foot and smooth step rollover.
- A unique structural bond connecting the two blades that eliminates fasteners and maintains flexibility to the toes. It also allows the two blades to flex and move independent of each other, increasing the ability to deflect in multiple directions, keeping a stable base of support for the user on uneven surfaces.
- A high energy return rubber heel shock that absorbs vertical loads and provides a fulcrum point for the Half Rocker Sole Plate to deflect around to provide a smooth forward progression over the foot.
- A technical fiber sock which limits noise and wear between the foot shell and the composite foot.
- A foot shell, available in dark or light color, which transmits the mechanical stresses between the shoe and the composite structure.
- A wedge to adjust heel stiffness.

2—Contents at Arrival

1. RUSH Foot[®]
2. Foot Shell
3. Technical Fiber Sock
4. Heel Wedge

3—Indications for Use and Maximum Weight

The HiPro, RAMPAGE, and RAMPAGE LP feet are designed exclusively for lower limb prostheses.

They are appropriate for patients with low, moderate, or high impact activity levels whose weight (carried load included) complies with the table below.

Advice: Patients of similar weight and size often require different category selections due to the activities they perform in their day to day life. Patients who perform high impact activities such as sprinting and basketball benefit from a slightly stiffer and more responsive foot compared to a patient who only performs low activities such as walking and playing golf. Proteor USA has created four different activity levels with examples to aid Prosthetists in selecting the correct activity level for their patient.

The four activity levels and their examples are the following:

- Low: Walking, Golfing
- Moderate: Hiking, Skiing
- High: Baseball, Wakeboarding
- High Impact: Sprinting, Basketball

The authorized maximum load is 365 lbs / 166 kg.

Each RUSH foot® is resistant to fresh, sea, and chlorinated water provided that the recommendations of Sections 8 and 9 are respected.

The HiPro, RAMPAGE, and RAMPAGE LP feet are designed for a 10mm heel height.

They include a male pyramid that allows the use of adult pyramid receivers.

To select the correct category for the patient, the Prosthetist must know the patient's activity level and weight. Using this information, the patient's category number can be found at the intersection of the patient's weight (vertical column) and the patient's activity level (horizontal column).

HiPro, & RAMPAGE / Rampage LP (x)										
Patient Weight	lbs	0 to 105	106 to 140	141 to 175	176 to 210	211 to 245	246 to 280	281 to 315	316 to 350	351 to 365
	kg	0 to 48	49 to 64	65 to 79	80 to 95	96 to 111	112 to 127	128 to 143	144 to 159	160 to 166
Low		1	2	3	4	5	6	7	8	9
Moderate		1	2	3	4	5	6	7	8	9
High		1 (2)	2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)
High Impact		2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)	- (-)

NOTE: Increasing category numbers indicate increased keel response.

FOR BILATERAL PATIENTS: Increase by one category for added stability

365lbs/166kg

4

3

2


1

4—Bench Alignment

Sagittal Plane:

After determining the appropriate socket flexion and heel height, position the weight line so that it falls through the center of the foot adapter. As the unique glass composite material provides substantially more flexibility than the other prosthetic feet, this represents a good starting point for a bench alignment.

Coronal Plane:

After determining the appropriate socket adduction/abduction, position the weight line so that it falls through the midline of the foot in a neutral M-L position.

5—Alignment (With the Patient)

Static Alignment:

Due to the shape of the rocker sole (bottom blade), patients may experience a new sensation when seeking the midpoint of the foot. The rocker sole allows patients to find their own comfortable static or standing position. The rocker bottom shape of the sole serves 2 primary functions :

- Provide a continuous and progressive point of contact throughout the entire step.
- Eliminate any "flat" or "dead" spot.

The Prosthetist is highly discouraged from introducing wedges at this stage of alignment.

Adjustment of the A/P set screws at the proximal adapter is the more appropriate place to make changes regarding plantar flexion or dorsiflexion, rather than the introduction of a wedge.

Dynamic Alignment:

The use of a slide adapter with the RUSH Foot® is highly encouraged for dynamic alignment purposes, as it provides the best solution to troubleshooting the most common alignment issues listed below:

- Hard or soft heel
- Hard or soft toe
- Varus or Valgus movements during Stance Phase.

Once the optimal relative socket/foot alignment is determined, RUSH Foot® encourages the use of plantar flexion / dorsiflexion at the proximal foot adapter to optimize heel-to-toe comfort and energy return.

The patient should wear the RUSH Foot® for at least one week to fully experience the responsiveness of a non-carbon fiber device.

If the patient still requests additional heel stiffness, it is appropriate to introduce and utilize the wedge provided in the original RUSH Foot® package.

This wedge uses a double-sided sticker so the lower face of the blade should be degreased before using. See figures on the next page for correct wedge placement.

Proper Wedge Placement

Improper Wedge Placement

Placement between components will adversely affect foot performance and void warranty.

6—Finishing

To avoid rubbing noises, pull the sock up to the base of the pyramid and immobilize it using a fine string for example.

7—Product Life

This component has been tested cyclically, which corresponds to at least 3 years of use depending on the activity of the patient.

An annual check by the Prosthetist is recommended.

8—Maintenance

No maintenance operation such as lubrication, work on the screws, or other parts is required for this foot.

The patient is encouraged to periodically remove the foot shell and clean the foot to remove any dirt or debris which could alter the foot behavior or damage it.

To clean the foot shell, use a slightly damp cloth or sponge.

Do not attempt to remove the foot from its shell by pulling it out by hand. This may damage the foot, voiding the warranty.

Always apply or remove the foot shell with an appropriate tool or device to minimize wear and damage, such as a shoe horn.

The foot shell and the sock are the only foot components that may need to be replaced (if worn out, lost, etc.).

If the foot shell and/or technical fiber sock need to be replaced, please contact Proteor USA for a replacement. The RUSH Foot® is designed to work exclusively with the Proteor USA foot shell and supplied technical fiber sock. Only the Prosthetist can order a replacement foot shell or technical fiber sock.

9—Advice for Use and Safety

A technical fiber sock should always be used with the RUSH Foot® to minimize or eliminate the introduction of debris that may come into contact with the composite material, as well as minimize wear to the foot shell.

Use temperature : - 20°C / + 60°C (- 4°F / 140°F) .

Use of the foot that does not adhere to your Prosthetists recommendations can damage foot components (carrying heavy loads for example).

In case of an abnormal behavior of your foot, or if you feel changes in your foot features, or if the foot sustained a hard impact, contact your Prosthetist.

After use in water or in a very dirty environment:

- Remove the foot shell and take off the technical fiber sock.
- Rinse the foot and the other parts thoroughly with clean water to eliminate debris such as dirt, sand, rock, etc.
- Allow all the components to dry before reassembly.

Failure to follow these guidelines may result in premature failure and void the foot warranty.

10—Recycling

The different items of the foot are special wastes and must be handled according to local laws.

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Proudly handcrafted in the U.S.A.
PROTEOR SAS