
OUTSIDER by HOPPER®

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EN

EN : OUTSIDER by HOPPER®

Intended for orthoprosthethists, this manual is designed to optimize the adjustment of the **Outsider by Hopper®** prosthetic foot and to assist in selecting the most suitable configuration for the user.

PRODUCT DESCRIPTION

The Outsider by Hopper® foot is a medical device intended for the prosthetic fitting of individuals with lower-limb amputations. It is designed to provide a versatile solution combining energy return, comfort, and stability, enabling users to engage in a wide range of physical activities while remaining suitable for an active daily lifestyle. The maximum authorized user weight is 130 kg.

Device weight (including adapter and modules) : **570 to 685 g**

Size range : **24 to 30 cm**

Suitable for users weighing **35 to 130 kg**

Suitable for a **wide range of activities** and **waterproof**

ASSEMBLY INSTRUCTIONS

- 1 Select the appropriate reference for the user


Weight [kg]	35-50	51-65	66-80	81-95	96-110	111-130
Weight [lbs]	77-111	112-143	144-176	177-209	210-243	244-287
Size 24-25 [cm]	TXS-A	TS-A	TM-A	TL-A	TXL-A	
Size 26-30 [cm]	TXS-B	TS-B	TM-B	TL-B	TXL-B	TXXL-B
Footshell Modules	Sh1 <95kg	Sh1 <95kg	Sh1 <95kg	Sh1 <95kg	Sh2 >95kg	Sh2 >95kg

- 2 The heel module (A) must be cut along the lines indicated on the bottom side of the system, according to the desired size (e.g. for size 27 below).


⚠ NEVER CUT THE COMPOSITE BLADE


- 2bis For sizes 24 and 25 only, the forefoot module (B) can be cut along the marked line.


- 3 Using the hook-and-loop fastener, assemble the modules with the blade


- 4 Add a sock to make fitting easier and protect the modules

- 5 Put on the shoe and tighten the laces **as much as possible**, starting from the bottom

6


The alignment reference line must:

- Pass through the center of the socket (X)
- Reach the center of the adapter (Y), keeping the tube straight

7

For more dynamic use:

- Add plantar flexion.

⚡ SEE THE "DYNAMIC ADJUSTMENTS" SECTION FOR FURTHER GUIDANCE

PREPARING THE PRODUCT FOR USE

The installation of the product, as well as socket adjustments, must be carried out by a qualified healthcare professional and performed in accordance with the instructions relating to alignment, assembly, and adjustment procedures.

CHOICE OF DEVICE

Weight [kg]	35-50	51-65	66-80	81-95	96-110	111-130
Weight [lbs]	77-111	112-143	144-176	177-209	210-243	244-287
Size 24-25 [cm]	TXS-A	TS-A	TM-A	TL-A	TXL-A	
Size 26-30 [cm]	TXS-B	TS-B	TM-B	TL-B	TXL-B	TXXL-B
Footshell Modules	Sh1 <95kg	Sh1 <95kg	Sh1 <95kg	Sh1 <95kg	Sh2 >95kg	Sh2 >95kg

WARNING

- The Outsider foot must be used in combination with Outsider footshell modules.
- Exceeding the weight limit may result in system failure.

ASSEMBLY INSTRUCTIONS

- [see Fig. 2 & 2bis] Trim the modules to dimensions appropriate for the user. We recommend first making a straight cut, followed by the rounded cuts. If the modules were supplied pre-cut to the patient's shoe size, skip this step. **WARNING:** Never cut the composite blade.
- [see Fig. 3] Install the heel module (A) under the blade by positioning the end of the module (the side where the size is indicated) flush with the side edge of the blade. Then add the forefoot module (B).
- [see Fig. 4] Place a sock around the complete system before inserting it into the shoe, in order to facilitate insertion and protect the integrity of the modules.
- [see Fig. 5] Tighten the laces as firmly as possible, starting from the bottom eyelets upward (to ensure even lacing across the entire area and create sufficient compression to secure the modules firmly inside the shoe).

WARNING

- The use of overly soft shoes with insufficient support is not recommended. This may result in reduced support and stability, as well as premature wear of the modules.
- In the case of combination with other prosthetic components, ensure, by referring to the respective instructions for use of each prosthetic element, that their compatibility and combined use are approved.

INFORMATION : TRANSFEMORAL AMPUTATION

- For transfemoral assembly with a prosthetic knee, refer to the manufacturer's instructions for the knee being used.
- When assembling the Outsider prosthetic foot with a prosthetic knee, it is recommended to prioritize the knee manufacturer's alignment and assembly guidelines. Indeed, the dynamics of our foot are oriented forward and may influence the behavior of the knee. This may affect the timing of the unlocking mechanism and could surprise the user. Ensure that the user is informed and safely supported during their first steps.

HEIGHT ADJUSTMENT BEFORE STATIC FITTING

The foot's design involves vertical compression of the blade to provide comfort and dynamism for the user. To compensate for these effects, add 5 to 10 mm to the socket-to-ground measurement.

STATIC ADJUSTMENTS

The main objectives of the static adjustment phase are:

- Ensure even weight distribution when the patient is standing on both legs.
- Check prosthetic height, internal/external rotation, and stability (anterior/posterior).

ALIGNMENT INSTRUCTIONS	
STEP 1	Check the appropriate length of the prosthesis [see height adjustment instructions]
STEP 2	Heel height: once the shoe is installed on the foot, check that the prosthetic tube is vertical (for walking use).
STEP 3	External rotation of the foot: approx. 5° Take into account the individual abduction or adduction position of the residual limb.
STEP 4	Sagittal orientation: follow the alignment line of the prosthetic foot [see fig. 6]
STEP 5	Connect the prosthetic foot and the prosthetic socket using the selected adapters. During adjustment and assembly, follow the instructions for use provided for the adapters.
STEP 6 [see fig. 6]	• Determine the centre of the socket using the 50/50 template. Position point X (socket midpoint) • Take into account the patient's individual flexion. • Ensure that the vertical alignment between point X (socket midpoint) and point Y (pyramid centre) is maintained.

EN

DYNAMIC ADJUSTMENTS

Adjustments during use are recommended in order to best adapt the settings to the specific needs and requirements of the user. Before making any adjustments, ensure that the patient understands how the system works.

The smoothness of the transition from heel to forefoot and comfort can be influenced by the following parameters:

- The stiffness of the heel or carbon blade
- The anterior/posterior position of the foot
- Plantar flexion or dorsiflexion of the foot
- The type of shoe

If necessary, consider the following actions:

SYMPTOMS	ACTIONS
The patient may feel as if they are collapsing at the end of the step [see fig. 7].	Check that the theoretical alignment is respected [see fig. 6]. The anterior translation of the socket may be too significant. Consider increasing plantar flexion. Consider using a stiffer blade.
Extra effort is required to roll onto the forefoot and/or the forefoot is perceived as too stiff.	Check that the theoretical alignment is respected [see fig. 6]. The posterior translation of the socket may be too significant. Consider increasing dorsiflexion. Consider using a more flexible blade.
The foot is perceived as too "rigid." As an orthoprosthesis, visually, no deformation (compression) of the foot is observed. It is likely that the rollover is jerky and that the patient has to exert extra effort to complete the step.	The blade is too stiff, replace it with a more flexible blade.
The foot is perceived as too "soft." As an orthoprosthesis, you visually observe a vertical compression greater than 30 mm during walking, and the patient appears to collapse at the end of the step.	The blade is too flexible, replace it with a stiffer blade.

MAINTENANCE AND SAFETY

ENVIRONNEMENTAL CONDITIONS

- Humidity : The device is classified as aqua-compatible and can withstand contact with fresh water, chlorinated water, salt water, and sweat.

WARNING

- In case of immersion: Wearing a shoe is recommended to help keep the modules securely in place.

Occasionally permitted environmental conditions :

- Dust, sand, and dirt: refer to the cleaning instructions for proper maintenance.

USE WITHOUT SHOES

- The use of the device without shoes is not recommended, as it may lead to premature wear of the modules.
- Use of the device without shoes is strictly prohibited outdoors and on any rough surface

RISKS OF INJURY AND PRODUCT DAMAGE

- It is strictly forbidden to cut, drill, heat, or modify the product in any way.
- Discontinue use in the event of unusual pain, persistent discomfort, or loss of control, and consult a healthcare professional.
- Inspect the product regularly to ensure there are no signs of damage (wear, delamination).
- In the event of a malfunction or any doubt regarding the condition of the product, the patient must immediately stop using it and inform the manufacturer.

CLEANING

- Clean with fresh water after exposure to chemicals, sand, dust, or dirt, then dry using a soft cloth.
- Rinse with water and remove any sand, dust, or dirt from the foot and modules. Dry with a cloth after cleaning.
- After contact with salt or chlorinated water, rinse with clean water and dry using a soft cloth.

MAINTENANCE

- Schedule regular follow-up appointments with the patient, based on product use, to check settings and adjust them if necessary.
- Replace the modules when heel compression exceeds 5 mm.

EN

SERVICE LIFE

- **Composite blade**

No exact service life can be defined for use of the product during sporting activities, as it is subjected to highly variable stresses depending on the type and intensity of activity. The device has been tested for two million load cycles. Depending on the patient's activity level, this may correspond to approximately 2 to 3 years of use.

- **Foot Shell modules**

The modules are wear parts subject to normal wear and tear. No precise service life can be defined for use of the product during sporting activities, as the level and type of stress vary significantly depending on the activity and its intensity.

- **Hook-and-loop fastening system**

The hook-and-loop fastening system is subject to gradual wear directly related to the number of detachments.

WARNING: It is important to limit the number of times the modules are detached in order to preserve the effectiveness of the fastening system.

END OF LIFE OF THE DEVICE

It is recommended that the device be returned to the manufacturer at the end of its life so that it can be recycled or reused.

LIABILITY

The manufacturer recommends using the device only under the specified conditions and for the intended purpose. The device must be maintained in accordance with the instructions for use.

The manufacturer cannot be held liable for damage caused by combinations of components not authorized by the manufacturer.

COMPLIANCE

The Outsider by Hopper® feet comply with European Conformity (CE) certification standards. The composite blade and its upper adaptation system have undergone fracture testing and two million load cycle tests under conditions designed to closely replicate intended use.

Any serious incident occurring in relation to the device must be reported to Hopper SAS, as well as to the competent authority of the country in which the user is based.

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Dispositivo médico

Designation	Reference	UDI-DI
Outsider by Hopper TXS-A	HP01TXS-A	03770027689263
Outsider by Hopper TS-A	HP01TS-A	03770027689256
Outsider by Hopper TM-A	HP01TM-A	03770027689249
Outsider by Hopper TL-A	HP01TL-A	03770027689232
Outsider by Hopper TXL-A	HP01TXL-A	03770027689270
Outsider by Hopper TXS-B	HP01TXS-B	03770027689225
Outsider by Hopper TS-B	HP01TS-B	03770027689218
Outsider by Hopper TM-B	HP01TM-B	03770027689201
Outsider by Hopper TL-B	HP01TL-B	03770027689195
Outsider by Hopper TXL-B	HP01TXL-B	03770027689188
Outsider by Hopper TXXL-B	HP01TXXL-B	03770027689171
Outsider Footshell multiaxial upper and heel pads Sh1 <95kg Size 24	HM02-1-24	03770027689430
Outsider Footshell multiaxial upper and heel pads Sh1 <95kg Size 25	HM02-1-25	03770027689423
Outsider Footshell multiaxial upper and heel pads Sh1 <95kg Size 26	HM02-1-26	03770027689416
Outsider Footshell multiaxial upper and heel pads Sh1 <95kg Size 27	HM02-1-27	03770027689409
Outsider Footshell multiaxial upper and heel pads Sh1 <95kg Size 28	HM02-1-28	03770027689393
Outsider Footshell multiaxial upper and heel pads Sh1 <95kg Size 29	HM02-1-29	03770027689386
Outsider Footshell multiaxial upper and heel pads Sh1 <95kg Size 30	HM02-1-30	03770027689379
Outsider Footshell multiaxial upper and heel pads Sh2 >95kg Size 24	HM02-2-24	03770027689362
Outsider Footshell multiaxial upper and heel pads Sh2 >95kg Size 25	HM02-2-25	03770027689355
Outsider Footshell multiaxial upper and heel pads Sh2 >95kg Size 26	HM02-2-26	03770027689348
Outsider Footshell multiaxial upper and heel pads Sh2 >95kg Size 27	HM02-2-27	03770027689331
Outsider Footshell multiaxial upper and heel pads Sh2 >95kg Size 28	HM02-2-28	03770027689324
Outsider Footshell multiaxial upper and heel pads Sh2 >95kg Size 29	HM02-2-29	03770027689317
Outsider Footshell multiaxial upper and heel pads Sh2 >95kg Size 30	HM02-2-30	03770027689300

To find more content about the Outsider foot and instructions for use in other languages, scan the following QR code: