

Exalto

ENTHUSIASM EVERY DAY



Carbon fiber energy storing foot

Designed for a daily safe and comfortable usage of patients with low to moderate activity.



20222212-UK-PL-00106-EXALTO

#HUMANFIRST

Exalto

ENTHUSIASM EVERY DAY



With its light weight optimized shape, and low built height, Exalto suites a large patient panel, up to 166 kg.



Two split blades, both covering the full foot length, to ease stability and comfortable rollover.



Resistant to immersion in salted and chlorine water.

Adapts to uneven terrains thanks to its split blades enabling inversion/eversion movements.





Exalto

ENTHUSIASM EVERY DAY



F13-00-_ A_ -RU

Category
01-09

Size
22-31

Size (cm)	Heel height	Build height	Weight
22 to 25	10 mm 3/8"	98 mm / 3 7/8"	440 g / 1 lbs
26 to 28		110 mm / 4 3/8"	540 g / 1.2 lbs
29 to 31		118 mm / 4 3/4"	700 g / 1.5 lbs

Selection of foot category based on patient's weight and activity										
Weight lbs	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
Weight kg	44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
Low	01	01	02	03	04	05	06	07	08	09
Moderate	01	02	03	04	05	06	07	08	09	-

Based on experience, it is recommended to increase the resistance for patients whose weight/activity is at the limit of the correspondence table.

FTC-3M-O_ _4R_

Size
22-31

Color
L- Light
M- Medium
D- Dark

Side
L- left
R- right



SPECTRA SOCK	Item n°.
22 à 25	SO-NPS-20030-00
26 à 28	SO-NPS-20035-00
29 à 31	SO-NPS-20040-00



BUMPER	Item n°.
	KIT-00-1148U-0