



RISE AND RUN.



In the gym, at the park, or on the trail... the EASY RUN blade by Proteor USA is the perfect foot to get you moving.

So, you're not an experienced sprinter? Not to worry! Whether you're exploring new activities or looking to take your daily workout to a new level, the **EASY RUN** will help you increase strength, mobility, happiness, and amplify your life. Specifically designed for the "every day user", the **EASY RUN** is here to help you **#riseandEASYRUN**.

With its all-terrain, all-weather, curved fiberglass blade design, and anti-slip Vibram® Sole... you can take your **EASY RUN** workout outside with confidence.

EASY RUN is available in 6 vibrant designs... allowing you to choose one as unique as your personality! Replacement soles are available in both black and lime green to keep your **EASY RUN** on the move.

THINK STRONG. BE STRONG. FINISH STRONG.

Technical FEATURES

| | |
|--------------------------|---------------------|
| Weight | 2.40lbs / 1088 g |
| Build Height | 8.66in / 220 mm |
| Materials | Fiberglass + rubber |
| Rotatable Pyramid | Pyramid |

PATIENT WEIGHT AND ACTIVITY

| Weight | 88-108lbs 40-49kg | 110-130lbs 50-59kg | 132-152lbs 60-69kg | 154-174lbs 70-79kg | 176-194lbs 80-88kg | 198-218lbs 90-99kg | 220-242lbs 100-110kg |
|------------------------|----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-------------------------|
| Moderate Impact | | 1 | 2 | 3 | 4 | 5 | 6 |
| High Impact | 1 | 2 | 3 | 4 | 5 | 6 | |



**AVAILABLE
IN 6 VIBRANT
DESIGNS**

#riseandEASYRUN



1236 West Southern Avenue #101, Tempe, AZ 85282 • **855.450.7300**

proteorusa.com [f](#) [@](#) [v](#) [in](#) **#HumanFirst**

Human First

Human First





ROTATABLE PYRAMID ADAPTER

CURVED FIBERGLASS BLADE

LOWER BLADE WITH
OPTIMIZED CURVATURE

ANTI-SLIP SOLE



1236 West Southern Avenue #101, Tempe, AZ 85282 • 855.450.7300
proteorusa.com     #HumanFirst

Human First