

RUSH FOOT[®] EVAQ8 COLLECTION

Joel Goltry, BK, RUSH FOOT HiPro EVAQ8

**SIMPLISTIC.
EFFECTIVE.
EFFICIENT.**

Finally, an unprecedented elevated vacuum system that seamlessly integrates effective vacuum with the same rugged simplicity of all the **RUSH Foot[®]** designs. The **RUSH Foot EVAQ8 Collection** is the easy solution to providing maximum vacuum efficiency in just a few quick steps. No batteries or complicated, mechanical pump required.

Unlike the competition, upwards of 22inHg can be achieved with this surprisingly simple, integrated design that stabilizes volume, increases linkage, and heightens proprioception through elevated vacuum suspension. The **RUSH Foot EVAQ8 foot Collection** is all-terrain, fully-submersible (fresh or salt water), free of bulky components, and requires minimal maintenance. The **RUSH FOOT EVAQ8 Collection** is available in the **RAMPAGE, HiPro, ROGUE, and RAMPAGE LP** models providing the perfect option for the majority of individuals seeking elevated vacuum suspension.

**Go anywhere. Do anything.
The RUSH Foot EVAQ8
Collection elevates your life.**

RUSH ROGUE[®] EVAQ8 PROFILE INFORMATION

Weight Limit	365lbs./165kg	Heel Rise	3/8"/10mm	Size	Build Height*
Size	22cm - 29cm	Warranty	36 Months	22 - 24cm	6.5"
Foot Shell	Light or Dark	†Suggested Codes	L5987, L5984 L5781	25 - 27cm	6.8"
Product Weight	1.93lbs./875g (26cm)			28 - 29cm	7.3"

RUSH HiPro[®] EVAQ8 PROFILE INFORMATION

Weight Limit	365lbs./165kg	Heel Rise	3/8"/10mm	Size	Build Height*
Size	22cm - 29cm	Warranty	36 Months	22 - 24cm	6.5"
Foot Shell	Light or Dark	†Suggested Codes	L5981, L5781	25 - 27cm	6.8"
Product Weight	1.70lbs./770g (26cm)			28 - 29cm	7.3"

RUSH RAMPAGE LP EVAQ8 PROFILE INFORMATION

Weight Limit	365lbs./165kg	Heel Rise	3/8"/10mm	Size	Build Height*
Size	22cm - 29cm	Warranty	36 Months	22 - 24cm	4.5"
Foot Shell	Light or Dark	†Suggested Codes	L5981, L5986 L5781	25 - 27cm	4.9"
Product Weight	1.74lbs. /788g (27cm)			28 - 29cm	5.0"

RUSH RAMPAGE EVAQ8 PROFILE INFORMATION

Weight Limit	365lbs./165kg	Heel Rise	3/8"/10mm	Size	Build Height*
Size	22cm - 29cm	Warranty	36 Months	22 - 24cm	6.5"
Foot Shell	Light or Dark	†Suggested Codes	L5987, L5986 L5781	25 - 27cm	6.9"
Product Weight	1.31lbs/596g (26cm)			28 - 29cm	7.3"

† The listing of these codes should not be considered a guarantee for coverage or a payment. *Weights include footshell. * Build heights include footshell + 10mm heel rise.

CATEGORY CHART

Weight	lbs kg	0 - 105 0 - 48	106 - 140 49 - 64	141 - 175 65 - 79	176 - 210 80 - 95	211 - 245 96 - 111	246 - 280 112 - 127	281 - 315 128 - 143	316 - 350 144 - 159	351 - 365 160 - 166
Low Activity	1	2	3	4	5	6	7	8	9	
Moderate Activity	1	2	3	4	5	6	7	8	9	
High Activity	1 (2)	2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)	
High Impact	2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)	- (-)	

NOTE: Increasing category numbers indicate increased keel response. FOR BILATERAL PATIENTS: Increase by one category for added stability. (RAMPAGE LP)

Call us toll-free at 855.450.7300 or visit us at proteorusa.com

f t YouTube in p #HumanFirst



HumanFirst