



INTRODUCING THE GERY...

a comfortable, re-engineered low K2 foot.



The newly re-designed GERY foot aims to provide lower level K2 users with the freedom they deserve.

- Greatly assists your patient's ability to ambulate with it's lightweight design.
- Accommodates heel heights ranging from .2" to .8" (5mm to 20mm).
- Cleaned easily with water and is resistant to dirt thanks to its closed cells EVA foam.



GET BACK TO YOUR LIFE.

- Engineered specifically for the low level K2, AK & BK population. Rated up to 150kg / 330lbs.
- Re-engineered with a new core
- 1 Year warranty
- 100 % Waterproof
- Comfort thanks to its design that reduces residual limb pain. Lightweight thanks to its EVA structure \approx 30% lighter than a SACH foot! Easy step completion thanks to its profiled core.
- Comes with malleolar clip

TECHNICAL FEATURES

Order number	1A200H	1A200BH
Weight	.87lbs. (395g)	
Build Height	3.22" (82mm)	
Structure material	EVA	

Heel height (T)	from .2" to .8" (5mm to 20 mm)
Connection	Pyramid
Maximum patient's weight	330lbs. (150kg)
Activity level	1-2

REF.	Side	Size	Suggested Codes	Color
1A200-H	G / D Left / Right	22cm to 29cm	L5972 - L5986	Caucasian
1A200-BH				Ebony

